




RAFTER D GENETICS

INTERNATIONAL EMBRYO TRANSFER CONSULTING

Embryo Transfer Training

Embryo Transfer Service

Semen & Embryo Sales




Contact Dr. Clifford Dorn
Rafter D Genetics
7750 Raymond Stotzer Pkwy
College Station, TX 77845
(979)260-7852
cgdorn@gte.net

Beltron EFT3002 Embryo Freezer

The EFT-3002 Embryo Freezer sets on top of your existing Semen tank and uses your LN2 to provide the cooling. Once the lower rod is frozen, there is no additional consumption of LN2 from your tank. The freezer can operate using 110v AC or 12v DC.



The freezing chamber will accommodate 64 straws in 32 drilled chambers. The freezing rate can be adjusted from 0.1°C/min to 0.6°C/min. The starting temperature is preset at -6.5°C and the final temp is -34.9°C. Comes with a carrying case & 12V adapter.



\$ 3,800.00 US

We also have used Nikon stereoscopes, new compound scopes & ultrasounds.

Contact Rafter D Genetics for info & Ordering!

(979)260-7852

<http://rafterdgenetics.com>